




Miltenyi Biotec

Changing the way of generating engineered T cells

CliniMACS Prodigy® T Cell Transduction Process



CliniMACS Prodigy® T Cell Transduction Process

It is happening – a new era of cell therapy

Chimeric antigen receptor (CAR) T cell therapy is undoubtedly leading a revolution in cancer therapy^{1,2}. The recent successes of CAR T cell therapy in fighting hematologic malignancies have led to a tremendous increase in interest in the immunotherapeutic field³. The great potential of genetically modified T cells now even expands into the area of solid tumors and infectious diseases⁴.

Currently, the manufacturing process of engineered T cells consists of various complex procedures, is labor intensive, and represents one of the biggest challenges in this area.

The CliniMACS Prodigy® T Cell Transduction (TCT) Process provides a unique all-in-one solution to this challenge. The CliniMACS Prodigy Platform enables the generation of gene-modified T cells in a standardized and fully automated way. It is a breakthrough in cell manufacturing, a leap to the next generation of engineered cell manufacturing^{5,6}.

The CliniMACS Prodigy TCT Process, a complete manufacturing solution – ready for the future.



References

1. Anurathapan, U. *et al.* (2014) *Cytotherapy* 16: 713–733.
2. Maus, M.V. *et al.* (2014) *Blood* 123: 2625–2635.
3. Barrett, D.M. *et al.* (2015) *J. Immunol.* 195: 755–761.
4. Whilding, L.M. and Maher, J. (2015) *Mol. Oncol.* 9: 1994–2018.
5. Mock, U. *et al.* (2016) *Cytotherapy* 18: 1002–1011.
6. Priesner, C. *et al.* (2016) *Hum. Gene Ther.* 27: 860–869.



Benefits of the TCT Process

Automation

- Simplifying a complex procedure.

Robustness

- Standardized procedure.

GMP compliance

- High-quality raw materials available.

Closed system

- Reliable product safety.
- Reduced hands-on time.

Transferable and scalable

- Supports centralized and de-centralized manufacturing models.

VIDEO



See the CliniMACS Prodigy® TCT Process in action!

Generate gene-modified T cells in a simple and automated fashion. Easy to use, this unique process will surely change the way you work.



► miltenyibiotec.com/tct

Automated workflow at a glance


Ready-to-use process to manufacture engineered T cells

The CliniMACS Prodigy® TCT Process enables the automated generation of engineered T cells. T cells are selected, transduced, expanded and the end product is reconstituted in its final formulation. Various starting materials can be used, such as leukapheresis or whole blood (WB). Final yields of up to 3×10^9 T cells are generated.

Specifications of the CliniMACS Prodigy® TCT Process

Starting material and volume		Final product	
Total cells	up to 20×10^9 total leukocytes	Total cells	on average 3×10^9
Volume	50–280 mL	Elution volume	100 mL
Material	WB, PBMCs, apheresis	Process time	10–14 days
Max. cell number for labeling:	3×10^9	Hands-on time	about 4 h (including material preparation)

Day	Step	Hands-on time at the instrument
Day 0	Cell enrichment	1 h
Day 0	Cell activation	20 min
Day 1 or 2	Cell transduction	10 min
Day 2–10	Cell expansion	10 min
Day 10–14	Final formulation	15 min



Cell enrichment

Select the best

To obtain robust, standardized, and reproducible results in the manufacture of engineered T cells, the key is to start with a well-defined enriched T cell population:

- **Efficiency:** Improves transduction performance
- **Cost-saving:** Avoids vector consumption by unwanted cells
- **Reproducibility:** Reduces variability of final cellular products between manufacturing runs

T cell subsets can be enriched by a clinical-scale magnetic cell separation system. Process-compatible reagents are:

- CliniMACS® CD4 Reagent
- CliniMACS CD8 Reagent
- CliniMACS CD62L Reagent




Figure 1: Flow cytometric analysis shows purity of CD4⁺ and CD8⁺ T cells enriched in an automated fashion by the CliniMACS Prodigy® TCT Process.




Figure 2: Composition of the cell populations after an automated co-enrichment with CliniMACS CD4 and CD8 Reagent on the CliniMACS Prodigy. Results from samples obtained from healthy donors (top) or patients (bottom) are shown. Each symbol represents one run on the CliniMACS Prodigy.

Cell activation

Effective activation and expansion

MACS® GMP T Cell TransAct™, with its unique features, enables biologically appropriate activation for your T cell manufacturing.

- **Unique format:** A colloidal polymeric nanomatrix conjugated to humanized, recombinant CD3 and CD28 agonists.
- **Ease of use:** Effective stimulation for a large range of cell densities.
- **Convenience:** Excess reagent is simply washed off.




Figure 3: CD4⁺ and CD8⁺ T cells from a healthy donor (upper row) or a melanoma patient (lower row) were enriched and cultured in TexMACS™ GMP Medium supplemented with IL-7 and IL-15 and activated using the TransAct™ T Cell Reagent. Pictures were taken with the integrated microscope camera of the CliniMACS Prodigy® 24 and 72 hours after activation.




Figure 4: The optimal transduction efficiency is obtained on day 1 after the addition of T Cell TransAct. Lentiviral vector encoding GFP was added to CD4⁺ and CD8⁺ T cells, co-enriched with CliniMACS® CD4 and CD8 Reagents, at different time points. Cells were cultured in TexMACS Medium supplemented with IL-7 and IL-15 and in the presence of T Cell TransAct. The transduction efficiency was determined on day 14 by flow cytometric analysis.

Cell transduction

Seamless integration

To achieve long-term gene expression through stable genomic insertion two types of gene transduction tools are commonly used: gamma-retroviral vectors (RV) or lentiviral vectors (LV). The CliniMACS Prodigy® TCT Process supports transduction with both vectors.

For easy integration of RV-based transduction into the automated workflow of the CliniMACS Prodigy TCT Process, Vectofusin-1® enhances your RV transduction efficiency. Vectofusin-1 is compatible with both static and spinoculation protocols.

For seamless integration into the automated CliniMACS Prodigy TCT Process we deliver tailored LV solutions. Customers profit from scalable LV manufacturing that meets US and EU GMP requirements and from a strong expertise in advancing LV design and technologies.




Figure 5: Transduction rates of enriched CD4⁺ and CD8⁺ T cells transduced on day 2 with gamma-retroviral GFP vector (GALV) could be increased with or without Vectofusin-1. Data shown for transduction performed in a small scale or on the CliniMACS Prodigy. This effect is further increased in an automated TCT Process on the CliniMACS Prodigy.




Figure 6: Enriched CD4⁺ and CD8⁺ T cells stimulated with MACS® GMP T Cell TransAct™ were transduced on day 1 with a LV encoding CD20 CAR. (A) Transduced T cells from healthy donor or lymphoma patients were analyzed by flow cytometry, (B) comparable transduction efficiencies of CD20 CAR-modified T cells were observed (healthy donors: n = 9; patients: n = 5).

Cell expansion

Culture is key

Optimal cultivation and expansion of transduced T cells rely on the strong synergy of MACS® GMP T Cell TransAct™ (page 6), TexMACS™ GMP Medium, and MACS GMP Cytokines.

TexMACS™ GMP Medium

- Specifically developed for T cell cultivation
- Serum- and xeno-component free
- Pharmaceutical-grade human serum albumin
- QC functionality test on every batch

MACS® GMP Cytokines

- Lot-to-lot consistency and lot-specific certificates of analysis
- Designed according to the recommendations of <USP 1043> on ancillary materials
- Manufactured and tested under a certified ISO 13485 quality system

Cell expansion takes place in the CentriCult Unit (CCU), an integrated cell cultivation chamber, which is part of the CliniMACS Prodigy® Tubing Set in order to maintain a closed system. The CCU allows culture volumes of 50 mL to 250 mL with an expansion capacity of up to 2×10^7 cells/mL. The culture duration is flexible and programmable, along with medium supply and exchange. The CCU also enables different agitation settings to maintain optimal culture conditions.




Figure 7: Enriched CD62L⁺ cells were cultured in the CliniMACS Prodigy during the TCT Process. Cell viability, pH-value, and glucose concentration throughout the whole process are shown.




Figure 8: Enriched CD4⁺ and CD8⁺ T cells were automatically expanded in the CliniMACS Prodigy after polyclonal stimulation, and transduction. Material from either healthy donors (n = 12) or patients (n = 6) was used. The T cell culture was monitored at different time points to determine the viability (A). The absolute cell count of viable T cells was calculated (B).

At the end of the CliniMACS Prodigy TCT Process, the final product is harvested in 100 mL buffer of your choice. Gene-engineered T cells are then ready for cryopreservation or direct transfusion.

- Up to 3×10^9 T cells
- Robust, reliable, and reproducible results
- High cell viability
- Automated harvest

In-process and quality control

Reliable reproducibility

In-process control (IPC) and quality control (QC) are required for reliable cell manufacturing. As an ideal partner for IPC/QC, Miltenyi Biotec with its MACS® Flow Cytometry portfolio provides complete solutions with unmatched convenience. Reach out to your local Miltenyi Biotec flow cytometry specialists to obtain detailed information on flow antibody panels relevant for CAR T cell manufacturing. When using the MACSQuant® Instrument, Express Modes allow you to fully automate and standardize your flow cytometry process. The MACSQuant Instrument and the MACS Flow Cytometry portfolio are for research use only.




Figure 9: Phenotypes of enriched CD4⁺ and CD8⁺ cells immediately after enrichment (left) or after stimulation, transduction, and expansion (final product, right). The frequency of naive (TN: D62L⁺CD45RO⁻), central memory (TCM: CD62L⁺CD45RO⁺), and effector memory (TEM: CD62L⁻CD45RO⁺) T cells in the enriched population vs. the final product, manufactured in the CliniMACS Prodigy®, was analyzed by flow cytometry. Healthy donors, n = 13 (A). Patients, n = 5 (B).




Figure 10: Percentages of CD4⁺ and CD8⁺ T cells in the enriched population vs. the final product manufactured in the CliniMACS Prodigy were determined. Healthy donors, n = 13 (A). Patients, n = 5 (B).



Figure 11: The MACSQuant Analyzer 10 (RUO) is a powerful benchtop flow cytometer for highly sensitive, multiparameter cell analysis.

Functionality of manufactured cells

Effectiveness that matters

Gene engineered T cells are assessed via cytokine secretion and their tumor cell killing.




Figure 12: Longitudinal monitoring of tumor burden, mock GFP vs. CD20 CAR therapeutic effect from day 3 post therapy onwards. Mock-GFP fails to have any therapeutic effect; all animals had to be sacrificed by day 17.




Figure 13: Automatically manufactured CD20 CAR T cells are fully functional in terms of cytokine secretion. CD20 CAR T cells, derived from healthy donor (A, n = 5) or patient (B, n = 4) samples, were cocultured with JeKo-1 as the target cell line (effector:target = 1:1), and cytokine secretion was analyzed using the MACSplex Cytokine 12 Kit, human.




Figure 14: Cytotoxicity assay demonstrates an efficient and specific tumor cell killing by automatically manufactured CART cells. CD20 CAR T cells, which were derived from either healthy (n = 8) or patient material (n = 4), were co-cultured with JeKo-1 as the target cell line at the indicated effector: target cell (E:T) ratios for 24 hours. Percentage of tumor cell killing was determined by flow cytometry.

Consumables and software

Essentials for success

CliniMACS Prodigy® TS 520

- One single-use tubing set for the entire workflow
- Sterile filters at access ports
- A truly closed system
- Integrated sampling pouches allowing IPC/QC anytime

CliniMACS Prodigy® TCT Process

- Tailored for the manufacture of engineered T cells
- Enrichment of different T cell subsets possible
- Compatible with both RV and LV transduction applications
- GMP-compliant cytokines and cell culture media are available to support the manufacturing process



Figure 15: Tubing Set installed on the CliniMACS Prodigy Platform.




Figure 16: The complete process (indicated in dark blue) to generate genetically engineered T cells in a closed system, using the Tubing Set TS 520 on the CliniMACS Prodigy Platform. Spinoculation is an optional activity for the TCT Process.



Miltenyi Biotec

Germany/Austria/ Switzerland

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

USA/Canada

Miltenyi Biotec Inc.
2303 Lindbergh Street
Auburn, CA 95602, USA
Phone 800 FOR MACS
Phone +1 530 888 8871
Fax +1 877 591 1060
macs@miltenyibiotec.com

Australia

Miltenyi Biotec
Australia Pty. Ltd.
Unit 16A, 2 Eden Park Drive
Macquarie Park NSW 2113
Australia
Phone +61 2 8877 7400
Fax +61 2 9889 5044
macs@miltenyibiotec.com.au

Benelux

Miltenyi Biotec B.V.
Schipholweg 68 H
2316 XE Leiden
The Netherlands
macs@miltenyibiotec.nl

Customer service The Netherlands

Phone 0800 4020120
Fax 0800 4020100

Customer service Belgium

Phone 0800 94016
Fax 0800 99626

Customer service Luxembourg

Phone 800 24971
Fax 800 24984

China

Miltenyi Biotec Technology &
Trading (Shanghai) Co., Ltd.
Room 2309
No. 319, Xianxia Road
Changning District
200051 Shanghai, P.R. China
Phone +86 21 62351005
Fax +86 21 62350953
macs@miltenyibiotec.com.cn

France

Miltenyi Biotec SAS
10 rue Mercœur
75011 Paris, France
Phone +33 1 56 98 16 16
Fax +33 1 56 98 16 17
macs@miltenyibiotec.fr

Italy

Miltenyi Biotec S.r.l.
Via Paolo Nanni Costa, 30
40133 Bologna
Italy
Phone +39 051 6 460 411
Fax +39 051 6 460 499
macs@miltenyibiotec.it

Japan

Miltenyi Biotec K.K.
Nittsu-Eitai Building 5F
16-10 Fuyuki, Koto-ku,
Tokyo 135-0041, Japan
Phone +81 3 5646 8910
Fax +81 3 5646 8911
macs@miltenyibiotec.jp

Nordics and Baltics

Miltenyi Biotec Norden AB
Scheelevägen 17
223 70 Lund
Sweden
macs@miltenyibiotec.se
Customer service Sweden
Phone 0200-111 800
Fax 046-280 72 99

Customer service Denmark

Phone 80 20 30 10
Fax +46 46 280 72 99

Customer service

**Norway, Finland, Iceland,
and Baltic countries**
Phone +46 46 280 72 80
Fax +46 46 280 72 99

Singapore

Miltenyi Biotec Asia Pacific Pte Ltd.
100 Beach Road
#28-06 to 28-08 Shaw Tower
Singapore 189702
Phone +65 6238 8183
Fax +65 6238 0302
macs@miltenyibiotec.com.sg

South Korea

Miltenyi Biotec Korea Co., Ltd
Arigi Bldg. 8F
562 Nonhyeon-ro
Gangnam-gu
Seoul 06136, South Korea
Phone +82 2 555 1988
Fax +82 2 555 8890
macs@miltenyibiotec.co.kr

Spain

Miltenyi Biotec S.L.
C/Luis Buñuel 2
Ciudad de la Imagen
28223 Pozuelo de Alarcón (Madrid)
Spain
Phone +34 91 512 12 90
Fax +34 91 512 12 91
macs@miltenyibiotec.es

United Kingdom

Miltenyi Biotec Ltd.
Almac House, Church Lane
Bisley, Surrey GU24 9DR, UK
Phone +44 1483 799 800
Fax +44 1483 799 811
macs@miltenyibiotec.co.uk

www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS GMP Products are for research use and *ex vivo* cell culture processing only, and are not intended for human *in vivo* applications. For regulatory status in the USA, please contact your local representative. MACS GMP Products are manufactured and tested under a quality system certified to ISO 13485 and are in compliance with relevant GMP guidelines. They are designed following the recommendations of USP <1043> on ancillary materials. The CliniMACS System components, including Reagents, Tubing Sets, Instruments, and PBS/EDTA Buffer, are designed, manufactured and tested under a quality system certified to ISO 13485.

In the EU, the CliniMACS System components are available as CE-marked medical devices for their respective intended use, unless otherwise stated.

The CliniMACS Reagents and Biotin Conjugates are intended for *in vitro* use only and are not designated for therapeutic use or direct infusion into patients.

The CliniMACS Reagents in combination with the CliniMACS System are intended to separate human cells. Miltenyi Biotec as the manufacturer of the CliniMACS System does not give any recommendations regarding the use of separated cells for therapeutic purposes and does not make any claims regarding a clinical benefit.

For the manufacturing and use of target cells in humans the national legislation and regulations - e.g. for the EU the Directive 2004/23/EC ("human tissues and cells"), or the Directive 2002/98/EC ("human blood and blood components") - must be followed. Thus, any clinical application of the target cells is exclusively within the responsibility of the user of a CliniMACS System.

In the US, the CliniMACS CD34 Reagent System, including the CliniMACS Plus Instrument, CliniMACS CD34 Reagent, CliniMACS Tubing Sets TS and LS, and the CliniMACS PBS/EDTA Buffer, is FDA approved; all other products of the CliniMACS Product Line are available for use only under an approved Investigational New Drug (IND) application or Investigational Device Exemption (IDE). CliniMACS MicroBeads are for research use only and not for human therapeutic or diagnostic use.

In the US, the CliniMACS Prodigy T Cell Transduction Process is available for use under an approved Investigational New Drug (IND). In the US, the MACSQuant Instrument and the Flow cytometry portfolio are available for RUO.

In the US, unless otherwise indicated, Lentigen products are for research and investigational use, only. Investigational use is allowed under approved IND-application, only. For EU: Please check with your regulatory authority for specific guidance on your investigational product. CliniMACS, CliniMACS Prodigy, MACS, the MACS logo, MACSQuant, TexMACS, and TransAct are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only. Copyright © 2017 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.