

Towards optimal neutrophil purification

Alina Bartholomäus¹, Martin Mengel¹, Daryl Grummit², Christian Peth¹, and Martin Büscher¹ ¹Miltenyi Biotec GmbH, Bergisch Gladbach, Germany, ²Owl biomedical Inc., Santa Barbara, CA, USA

Introduction

The challenge of studying neutrophils is their short lifetime and their fragility. Therefore, there is a high demand for a fast and convenient purification protocol of neutrophils from whole blood. Currently, density gradient centrifugation is still the "gold standard" for neutrophil isolation, although it is time consuming and elaborate, and the achieved purities and yields are relatively low. Besides the isolation of human neutrophils from whole blood: i) flow sorting immunomagnetic separation using beads, flow sorting offers a promising alternative. However, it is difficult to sort viable and nonactivated neutrophils with conventional droplet-based sorters. The MACSQuant[®] Tyto[®] is a microchip-based cell sorter, which requires


no sheath fluid and droplet formation, thus applying very low shear forces and no decompression. With these characteristics, the MACSQuant Tyto offers gentle cell sorting conditions and is a wellsuited candidate instrument for neutrophil isolation. Here we compared the sort performance of three techniques for with the MACSQuant Tyto, ii) magnetic separation with the MACSxpress[®] Neutrophil Isolation Kit, and iii) Percoll[®] density gradient centrifugation. We analyzed viability, functionality, and activation status of the cells.

Results

Neutrophil separation

Human neutrophils were separated from whole blood with the three different methods. In the MAC SQuant Tyto sorting strategy, neutrophils were stained with CD45-APC and CD16-PE-Vio[®] 770, and SSC^{hi}CD16⁺ cells were sorted. The MACSxpress Neutrophil Isolation Kit is based on the depletion of non-neutrophil leukocytes by a cocktail of magnetic beads coupled to specific antibodies and the simultaneous sedimentation of erythrocytes. These two methods were compared to a standard density

gradient centrifugation protocol, which is composed of two steps. First, Histopaque[®] 1119 medium is used to separate granulocytes from PBMCs and erythrocytes. A subsequent discontinuous Percoll gradient enables further purification of neutrophils. All isolated neutrophil fractions were analyzed with the MACSQuant Analyzer 10 (fig. 1; a trigger was set on CD45-APC and viable cells among CD45⁺ cells were gated. Neutrophils were identified as SSC^{hi}CD16⁺ cells).




Determination of the optimal sort performance

In order to compare the different purification methods, we developed a metric to evaluate the sort performances with regard to the number of sorted cells and cell purity, as well as the time it took to process the sample (fig. 2B). With this sort performance index (SPI) it is now possible to quickly compare distinct experimental setups, different cell sorting instruments, or even different sorting technologies¹.

With a standard Percoll density gradient centrifugation, the achieved purity was around 93%, but only 30% of the target cells could be recovered and the protocol took around 1–1.5 h to

complete (table 1). The MACSxpress Neutrophil Isolation Kit allowed a high cell purity and yield, based on a 20-minute protocol. Neutrophils that were sorted with the MACSQuant Tyto showed a comparable purity. The yield however was lower and the separation time was around 1.5 h. Therefore, MACSxpress Separation had the highest SPI, whereas the SPIs of density gradient centrifugation and MACSQuant Tyto were in a similar range. Nevertheless, the MACSQuant Tyto protocol is much easier to perform, smaller blood volumes are needed, and higher purities can be achieved.






	MACSQuant Tyto	MACSxpress	Percoll
Purity (%)	95	97	93
Yield (%)	31	75	25
Processed cells/min	1.31×10 ⁵	9.00×10 ⁵	4.03×10 ⁵
Speed index*	7	47	22
Purity index**	65	78	59
Sort performance index	34	66	35
* Speed index = (Mean processed cells/min)/(Max. processed cells/min) ** Purity index = $[2 \times \log (100 - mean purity)] \times 50$			

Table 1


Cell viability and activation status

positive controls. There was no up-regulation of CD11b (fig. 3B), Cell viability was assessed by using annexin V as a marker for apoptotic cells and PI as a marker for dead cells. After 4 h of CD66b (fig. 3D), and CD63 (fig. 3E) in neutrophils isolated incubation in medium at 37 °C and 5% CO₂, neutrophils isolated with MACSxpress Technology and only minimal activation in cells isolated with the MACSQuant Tyto or density with MACSxpress Technology showed the highest viability, followed by cells sorted with the MACSQuant Tyto. Neutrophils gradient compared to the negative controls. In addition, no that were isolated with the density gradient showed the highest down-regulation of CD62L could be detected in neutrophils percentages of apoptotic and dead cells (fig. 3A). In addition, isolated with MACSxpress Beads (fig. 3C), and only a slight decrease in the MFI could be detected in cells that were different activation markers were analyzed. Fresh whole blood was used as negative control and PMA-stimulated, heparinobtained using a density gradient or the MACSQuant Tyto. treated or EDTA-treated whole blood samples were used as



Cell functionality

The ability to phagocytose particles and to produce reactive which is cleaved into Rhodamine 123 in the presence of an oxygen species (ROS) are hallmarks of neutrophils. Thus, oxidative burst (fig. 4B, top). All isolated neutrophils showed separated neutrophils were incubated with pHrodo[™] Red *E. coli* a low oxidative burst after fMLP and a high oxidative burst after BioParticles[®] and the percentage of cells that were able to PMA stimulation. The enzymatic activity was comparable to the phagocytose these particles was analyzed (fig. 4A). Although the control (fig. 4B, bottom). The ability to migrate along an fMLP percentage of phagocytic cells among isolated neutrophils was gradient was assessed in a 2D chemotaxis assay by using slides slightly lower than in the control, neutrophils isolated with all from Ibidi[®]. Figure 4C shows the traces of cell migration towards three techniques clearly showed phagocytic activity. a 1 µM fMLP source after 30 minutes. Cell migration was recorded with an IN Cell Analyzer microscope with a 10× objective. Frames ROS generation was assessed by stimulating the cells with *E. coli*, fMLP, or PMA in the presence of the fluorogenic substrate DHR123, were captured every 30 s.



Conclusion

- The MACSxpress Neutrophil Isolation Kit shows the best sort performance with regard to cell purity and yield and sorting speed. Cells displayed a high viability and preserved functionality.
- Both magnetic cell separation and fluorescence-based flow sorting of neutrophils with the MACSQuant Tyto offer a more convenient and sensitive purification protocol than the current gold standard, i.e., density gradient centrifugation.

References

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, the MACS logo, MACSQuant, MACSxpress, Tyto, and Vio are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. Percoll is a trademark of GE Healthcare companies. All other trademarks mentioned in this document are the property of their respective owners and are used for identification ourposes only. Copyright © 2018 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.

and a subsequent sort on the MACSQuant Tyto could provide a fast method to isolate neutrophil subsets.

A combination of the MACSxpress Neutrophil Isolation Kit

A. Bartholomäus *et al.* (2015) Proposal for a Sort-Performance Index. Poster Cyto 2015, Glasgow, Scotland.